

# Extras




## Recommendation

OASE filter foams are manufactured exclusively for use in pond water that is suitable for fish. Consequently, they exceed the special requirements. Strict quality controls and specifications ensure consistent and high quality workmanship of the filter foams.

## Tip!

Only regular replacement of the filter foams assures the continuous performance of your filter system.

## HOW TO FIND THE RIGHT SPARE PART FOR YOUR FILTER

	Item no.	Description	Scope of delivery		Filter	Required quantity
<b>Underwater filters</b>						
	13705	Set Filtral 2500	●		Filtral UVC 2500	1 x
	35836	Set Filtral 5000	●		Filtral UVC 5000	
<b>Pressure filters</b>						
	15558	Set BioPress 4000	● 2x ●		BioPress Set 4000	1 x
	15564	Set BioPress 6000 – 12000	● 2x ●		BioPress Set 6000	
					BioPress Set 8000	
					BioPress Set 10000	
			2x ● 2x ●		BioPress Set 12000	
	56109	Set FiltoClear 3000	● 2x ●		FiltoClear 3000	
	56110	Set FiltoClear 6000	3x ● 3x ●		FiltoClear 6000	
	56111	Set FiltoClear 11000			FiltoClear 11000	
	51255	Set FiltoClear 12000	2x ● 3x ●		FiltoClear 12000	
	56884	Set FiltoClear 15000			FiltoClear 15000	
51258	Set FiltoClear 16000	3x ● 4x ●		FiltoClear 16000		
51290	Set FiltoClear 20000/30000	4x ● 5x ●		FiltoClear 20000 FiltoClear 30000		
<b>Flow-through filters</b>						
	50901	Cartridge Set FiltoMatic 3000 / 7000			FiltoMatic CWS 7000	1 x
	50904	Cartridge Set FiltoMatic 6000 – 25000	● 2x ●		FiltoMatic CWS 14000 / 25000	2 x
	40456	Cartridge Set BioTec Premium 80000	2x ● 2x ● 2x ●		BioTec Premium 80000	1 x



### OASE quality criteria

- Regular eco-toxicological examinations ensure a permanent, 100 % compatibility for plants and animals in the pond, also for the sensitive micro-organisms as well.
- The pore size is permanently defined and consequently ensures a constant and regular flow of water through the foam.
- OASE filter foams are characterised by incredible durability at a high level of dimensional stability. This stability is extremely important, particularly for the cleaning processes in which the filters are compressed via the cleaning handles.

The prerequisite for guarantee claims under each guarantee\*<sup>p. 310</sup> is the exclusive use of original OASE spare parts!

	Item no.	Description	Scope of delivery		Filter	Required quantity
<b>Flow-through filters</b>						
	35791	Foam red BioSmart 5000–16000	●		BioSmart 5000	● ●
	35792	Foam blue BioSmart 5000–16000	●		BioSmart 7000 / 8000 BioSmart 14000 / 16000	2x ● ● 3x ● ● 3x ● ●
	56677	Foam red BioSmart 18000–36000	●		BioSmart 18000, BioTec 5.1	● ● ● 2x ● ●
	56678	Foam blue BioSmart 18000–36000	●		BioSmart 20000 / 24000 / 30000 / 36000, BioTec 10.1	● 4x ● ● 4x ● ●
	56679	Foam green BioSmart 18000–36000	●			● 4x ● ● 4x ● ●
	42893	Set red /violet BioTec 12 / 40000 / 90000	● ●		BioTec ScreenMatic <sup>2</sup> 12 / 40000 BioTec ScreenMatic <sup>2</sup> 90000	3x ● ● 3x ● ● 6x ● ● 6x ● ●
	42895	Set blue BioTec 12 / 40000 / 90000	2x ●		BioTec ScreenMatic <sup>2</sup> 12 / 40000 BioTec ScreenMatic <sup>2</sup> 90000	2x ● ● 4x ● ●
	42894	Set red /violet BioTec 18 / 36 / 60000 / 140000	● ●		BioTec ScreenMatic <sup>2</sup> 18 / 60000 BioTec ScreenMatic <sup>2</sup> 36 / 140000	3x ● ● 3x ● ● 7x ● ● 7x ● ●
	42896	Set blue BioTec 18 / 36 / 60000 / 140000	2x ●		BioTec ScreenMatic <sup>2</sup> 18 / 60000 BioTec ScreenMatic <sup>2</sup> 36 / 140000	2x ● ● 4x ● ●
	54030	Foam red BioTec 5 / 10 / 30	●		BioTec 5 / 10	● ● ●
	54031	Foam blue BioTec 5 / 10 / 30	●		BioTec 30	10x ● ● 4x ● ●
<b>Module filters</b>						
	26983	Foam blue wide ProfiClear M3	●		ProfiClear filter foam module	2x ● ●
	26984	Foam blue narrow ProfiClear M3	●			● ●
	27295	Foam red ProfiClear M3	●			3x ● ●
	27296	Foam black narrow ProfiClear M5	●		ProfiClear phosphate binder module	● ●
	27297	Foam black wide ProfiClear M5	●			5x ● ●